



Minimum Essential Media
w/ Earle's Salts w/ NEAA w/o L-Glutamine

CAT N°: L0430

Storage conditions: +2°C to +8°C, protected from light

Shelf life: 24 months

Composition: Displayed on website; also available on request

Colour: Red, orange solution

pH: 7.3 ± 0.3

Osmolality: 290 mOsm/kg ± 10%

Endotoxin: < 1 EU/ml

Sterility tests:

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Cell Growth test: Medium tested for the ability to support L929 cell growth.

Other tests: Not applicable

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro for research or further manufacturing only and not for use as an Active Pharmaceutical Ingredient or food or animal feed.

Application:

Minimum Essential Medium (MEM) with Earle's Balanced Salts is a modification of Eagle's earlier medium Basal Medium Eagle (BME), containing higher concentrations of the essential nutrients. These media have demonstrated the ability to support a variety of normal and transformed cells in culture and contain Earle's Balanced Salts, which make them suitable for use in atmospheres charged with CO₂ gas. This medium contains also non-essential amino acids which allow an optimal cell growth.

Uses:

Supplements, such as antibiotics, should be added as sterile supplements to the medium.

Storage conditions and shelf-life of supplemented product will be affected by the nature of the supplements. Sterile serum should not be re-filtered before or after being added to sterile medium because growth promoting capacity may be reduced upon re-filtration.

Add 10 ml/l of L-Glutamine 200 mM (CAT N° : X0550) before using this medium.

BioWest also recommend the addition of 10% of Fetal Bovine Serum in the medium.



Technical data sheet Media Liquid

Ref : FT.L0430an
Page : 2/2

Version date : 13/06/2022

Signs of Deterioration:

Medium should be clear and free of particulate and flocculent material. Do not use if medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.

Remarks: Not applicable