

# Lymphosep




## Best choice for your cells:

Lymphosep is a standardized and high-quality gradient density solution, designed for a simple and rapid isolation of mononuclear cells from human whole blood. Each batch is tested for bacteria, fungi, yeast and endotoxin.

- ◆ Standardized and high quality medium
- ◆ Simple & rapid isolation of lymphocytes from whole blood
- ◆ Optimal density of **1.077 g/ml**
- ◆ Produced and tested according to the highest quality standards (ISO 13485 and ISO 9001)

## Recommendation:

- ◆ Use diluted and treated blood (with anticoagulant) for optimum yield.
- ◆ For best results, use blood drawn less than 2 hours before manipulation.
- ◆ Do not use with blood drawn 24 hours or more before the manipulation.
- ◆ Do not use as an API or direct inoculation into human or animals.



## Ordering information:

Cat n°	Product name
L0560-100	Lymphosep, 100ml
L0560-500	Lymphosep, 500ml

Please visit [www.biowest.net](http://www.biowest.net) for more information

[www.biowest.net](http://www.biowest.net)

Rue de la Caille

49340 Nuaille

Tel.: (+33) 241 464 242

Email: [biowest@biowest.net](mailto:biowest@biowest.net)