

Fetal Bovine Serum (FBS) Extracellular Vesicles Depleted



Description:

Biowest can offer you the best alternative to reduce time, effort and variability in your experiments with Biowest Fetal Bovine Serum, E.V. Depleted.

- ◆ ≥ 95 % depletion of exosome (guaranteed minimum)
- ◆ Original animal serum depletion system
- ◆ Cell growth promotion validated
- ◆ No significant variation in physico-chemical analyzes
- ◆ Ideal solution for exosome & cell communication research
- ◆ No need for long lasting ultra centrifugation protocols

Product informations:



Frozen



-20°C



60 months

Ordering informations:

FBS E.V. Depleted, South America Origin

S181M - 050

50ml

FBS E.V. Depleted, Europe Origin

S141M - 050

50ml

Other origins available on request

Please visit www.biowest.net for more information.

Applications:

E.V. Depleted serum is treated with our proprietary ultrafiltration method. This treatment depletes the microvesicles naturally present in the serum without adding any compound to the serum.

The depletion rate is measured by ELISA method.

www.biowest.net

Rue de la Caille

49340 Nuaillé

Tel.: (+33) 241 464 242

Email: biowest@biowest.net