



Technical data sheet
Media Liquid

Ref : FT.X0551an
Page : 1/2

Version date : 30/08/2022

Glutamine stable 100X, 200mM

CAT N°: X0551

Storage conditions: Room Temperature

Shelf life: 24 months

Composition: L-Alanyl-L-Glutamine 43.444g/l

Colour: Colourless

pH: 6,0 ± 1,0

Osmolality: 200 ± 30 mOsm/kg

Endotoxin: < 10 EU/ml

Sterility tests:

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Cell Growth test: Medium tested for the ability to support L929 cell growth or SP2/O-Ag14 cell growth.

Other tests: Not applicable

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (e.g.: under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro for research or further manufacturing only and not for use as an Active Pharmaceutical Ingredient or food or animal feed.

Application:

L-Glutamine is an essential amino acid required by all mammalian and insect cells for their culture. It is a crucial component of many cell culture media and serves as a major energy source for cells in culture.

Glutamine stable (L-Alanyl-L-Glutamine), dipeptide derivative of L-Glutamine prevents the intramolecular cyclization reaction associated with solutions of L-Glutamine.

The dipeptide derivative is metabolized within the cells to yield L-Glutamine and L-Alanine. This results in more consistent delivery of L-Glutamine to cells and avoid toxic build-up of ammonia in the cells. Glutamine stable is therefore especially dedicated for ammonia-sensitive cell lines.



Technical data sheet
Media Liquid

Ref : FT.X0551an
Page : 2/2

Version date : 30/08/2022

Glutamine stable allows the formulation of cell culture media containing L-Glutamine that may be stored at 4°C for extended periods. Solutions containing these derivatives can be even autoclaved without appreciable degradation of the product (loss of the product < 5% after 30 minutes at 121°C).

Uses:

Supplement cell culture medium with appropriate volume to achieve desired concentration (Use the same concentration of Glutamine stable in your medium as normally required for regular L-Glutamine).

Signs of Deterioration:

Glutamin should be clear of particulates and flocculent material.

Do not use if Glutamin is cloudy or contains precipitate.

Other evidence of deterioration may include degradation of physical or performance characteristics.

Remarks: Not applicable