

COMPOSITION

Trypsin

	X0915 Liquid mg/l	X0920 Liquid mg/l	L0940 Liquid mg/l	X0930 Liquid mg/l	P0940 Powder mg/l	L0941 Liquid mg/l	
Inorganic Salts	EDTA Disodium Salt Dihydrate	/	/	254,8	2214	2548	372
	Potassium Chloride	200	400	200	4000	/	200
	Potassium Phosphate Monobasic Anhydrous	200	60	200	/	/	200
	Sodium Bicarbonate	/	350	/	/	/	/
	Sodium Chloride	8000	8000	7950	8000	8040	7950
	Sodium Phosphate Dibasic Anhydrous	1150	47,88	1150	/	/	1150
O.C.*	D-Glucose Anhydrous	1000	1000	1000	10000	/	1000
	Phenol Red Sodium Salt	/	/	/	/	/	/
	Trypsin 1:250	25000	25000	500	5000	500	/
	Recombinant Trypsin 2500 USP	/	/	/	/	/	50

* Other Components