

Fetal Bovine Serum IgG Depleted



BOVINE



SOUTH AMERICA




Your cells will feel the difference:

Some experiments require low levels of specific antibody. Fetal Bovine Serum (FBS) has natural levels of IgG that are too high for some cell culture applications. Biowest can offer you a perfect alternative without changing the performance of FBS.

- ◆ Performance is equivalent to non-treated FBS
- ◆ Ultra low IgG levels < 5 µg/ml
- ◆ No modification on the cell growth promotion
- ◆ No modification on the cell viability
- ◆ Ideal solution for cell culture and protein purification applications where naturally occurring levels of IgG are too high

Applications:

- ◆ Cell culture: hybridoma and macrophage culture
- ◆ Monoclonal antibody production
- ◆ Immunology assays



FBS IgG depleted has been treated with proprietary chromatography method.

Ordering information:

Cat n°

S181I-100
S181I-500

Product name

Fetal Bovine Serum (FBS) South America, IgG Depleted, 100 ml
Fetal Bovine Serum (FBS) South America, IgG Depleted, 500 ml

Please visit www.biowest.net for more information

www.biowest.net
Rue de la Caille
49340 Nuaillé
Tel.: (+33) 241 464 242
Email: biowest@biowest.net