

Recombinant Trypsin



LOW GLUCOSE



Animal component free cell detachment:

Biowest's Recombinant Trypsin is a genetically engineered protein expressed in *Escherichia coli*.

As such it is totally animal free and ready-to-use for cell culture applications. Free from contaminating enzymes like chymotrypsin as well as protein inhibitors. It is highly stable with a high purity > 95 %.

- ◆ Animal Free Origin
- ◆ No risk of virus, BSE or other adventitious agents
- ◆ High purity
- ◆ High stability over 2 years
- ◆ Ready-to-use solution (1X) with Recombinant Trypsin and EDTA, without Calcium, Magnesium and Phenol Red

Applications:

- ◆ Cell culture

Ordering information:

Cat n°	Product name
L0941-100	Recombinant Trypsin, 100ml

Please visit www.biowest.net for more information

www.biowest.net

Rue de la Caille

49340 Nuaillé

Tel.: (+33) 241 464 242

Email: biowest@biowest.net