

COMPOSITION

MC Coy's 5A Medium Modified

L0210 P0390
Liquid Powder
mg/l mg/l

	L0210	P0390	
Amino Acids	Glycine	7,51	7,51
	L-Alanine	13,36	13,36
	L-Arginine Monohydrochloride	42,14	42,14
	L-Asparagine Monohydrate	45,03	45,03
	L-Aspartic Acid	19,97	19,97
	L-Cysteine	24,24	24,24
	L-Glutamic Acid	22,07	22,07
	L-Glutamine	219,15	219,15
	L-Histidine Monohydrochloride Monohydrate	20,96	20,96
	L-Hydroxy-L-Proline	19,67	19,67
	L-Isoleucine	39,36	39,36
	L-Leucine	39,36	39,36
	L-Lysine Monohydrochloride	36,54	36,54
	L-Methionine	14,92	14,92
	L-Phenylalanine	16,52	16,52
	L-Proline	17,27	17,27
L-Serine	26,28	26,28	
L-Threonine	17,87	17,87	
L-Tryptophan	3,06	3,06	
L-Tyrosine Disodium Salt Dihydrate	26,1	26,1	
L-Valine	17,57	17,57	
Inorganic Salts	Calcium Chloride Anhydrous	99,95	99,95
	Magnesium Sulfate Anhydrous	97,69	97,69
	Potassium Chloride	400	400
	Sodium Bicarbonate	2200	/
	Sodium Chloride	6460	6460
	Sodium Phosphate Monobasic Anhydrous	504	504
Vitamins	Ascorbic Acid	0,56	0,56
	Choline Chloride	5	5
	D-Biotin	0,2	0,2
	D-Ca Pantothenate	0,2	0,2
	Folic Acid	10	10
	Myo-Inositol	36	36
	Nicotinamide (Nicotinic acid amide)	0,5	0,5
	Nicotinic Acid	0,5	0,5
	P-Aminobenzoic Acid (PABA)	1	1
	Pyridoxal Hydrochloride	0,5	0,5
	Pyridoxine Hydrochloride	0,5	0,5
	Riboflavin	0,2	0,2
	Thiamine Hydrochloride	0,2	0,2
Vitamin B12	2	2	
O.C.*	D-Glucose Anhydrous	3000	3000
	L-Glutathione Reduced	0,5	0,5
	Peptone from soybean	600	600
	Phenol Red Sodium Salt	11	11

* Other Components