

COMPOSITION

Trypsin

	CAS number	L0909 Liquid mg/l	L0910 Liquid mg/l	L0930 Liquid mg/l	L0931 Liquid mg/l	L0932 Liquid mg/l	
Inorganic Salts	EDTA Disodium Salt Dihydrate	6381-92-6	/	/	221,4	380	200
	Potassium Chloride	7447-40-7	200	200	400	400	400
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	200	200	/	60	60
	Sodium Bicarbonate	144-55-8	/	/	580	350	350
	Sodium Chloride	7647-14-5	8000	8000	8000	8000	8000
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	1150	1150	/	48	48
O.C.*	D-Glucose Anhydrous	50-99-7	1000	1000	1000	1000	1000
	Phenol Red Sodium Salt	34487-61-1	2	/	2	10	10
	Trypsin 1:250	9002-07-7	2500	2500	500	2500	2500