

CytogenAmnio Medium



Innovative solution for your cell culture operations:

CytogenAmnio (AGM) is a ready-to-use medium, specifically formulated for primary human amniotic cells and chorionic villi culture. AGM is to be used for human karyotyping and other prenatal cytogenetic diagnostic.

The product supports a very efficient cell growth and attachment, resulting in early chromosome analysis and in high quality karyograms.

The product contains pre-tested FBS and gentamycin.

Purpose of Use:

CytogenAmnio Medium product is used to optimize the culture conditions for amniotic cells, providing a nutrient-rich environment that supports their growth and development, which is crucial for prenatal diagnostics and biomedical research.

Application:

Human karyotyping
Human *in vitro* prenatal diagnostic

Ordering information:

Cat n°	Product name
AGM-100	CytogenAmnio Medium, 100ml
	<i>Other packagings available on request.</i>

Please visit www.biowest.net for more information

www.biowest.net

Rue de la Caille
49340 Nuaille

Tel.: (+33) 241 464 242

Email: biowest@biowest.net