



# CytogenLympho Medium



## Innovative solution for your cell culture operations:

CytogenLympho Medium (LGM) is a ready to use medium that is specifically developed for the cultivation of primary human blood lymphocytes, which are used in karyotyping and other cytogenetic postnatal diagnostic procedures.

With the optimized PHA concentration CytogenLympho Medium (LGM) supports higher mitotic index, excellent chromosome morphology and banding pattern for highest quality chromosome analysis.

The product contains Penicillin/Streptomycin and Phytohaemagglutinine (PHA).

## Purpose of Use:

CytogenLympho serves the cultivation and growth of human lymphocytes and is to be used exclusively for *in vitro* diagnostic purposes on samples taken from humans.

## Application:

Cell culture  
Human cytogenetic diagnostic

## Ordering information:

Cat n°	Product name
LGM-100	CytogenLympho Medium, 100ml
	<i>Other packagings available on request.</i>

[www.biowest.net](http://www.biowest.net)

Rue de la Caille  
49340 Nuillé

Tel.: (+33) 241 464 242

Email: [biowest@biowest.net](mailto:biowest@biowest.net)

Please visit [www.biowest.net](http://www.biowest.net) for more information