

Fetal Bovine Serum (FBS) IgG Depleted






Description:

Some experiments require low levels of specific antibody. Fetal Bovine Serum (FBS) has natural levels of IgG that are too high for some cell culture applications. Biowest can offer you a perfect alternative without changing the performance of FBS.

- ◆ Performance is equivalent to non-treated FBS
- ◆ Ultra low IgG levels < 5 µg/ml
- ◆ No modification on the cell growth promotion
- ◆ No modification on the cell viability
- ◆ Ideal solution for cell culture and protein purification applications where naturally occurring levels of IgG are too high

Product informations:

 Frozen	 -20°C	 60 months
---	--	--


Ordering informations:

FBS IgG Depleted, South America Origin	
S181I - 100	100ml
S181I - 500	500ml

Other packagings available on request
Other origins available on request

Applications:

Cell culture: hybridoma and macrophage culture
Monoclonal antibody production
Immunology assays



FBS IgG depleted has been treated with proprietary chromatography method.

Please visit www.biowest.net for more information.

www.biowest.net
Rue de la Caille
49340 Nuaille
Tel.: (+33) 241 464 242
Email: biowest@biowest.net