

Dulbecco's Phosphate Buffered Saline (DPBS) w/o Calcium w/o Magnesium



Description:

DPBS buffer is a saline solution used to maintain cells in vitro. It is commonly employed during cell passaging to remove trypsin inhibitors by washing, or during cell counting when dilution is required.





DPBS is also used as a dilution buffer for reagents and/or samples in protein and cellular analysis assays, such as Western blotting, immunocytochemistry, and ELISA.

This formulation is calcium- and magnesium-free.

Applications:


- | | |
|--------------------|--------------------------------|
| Cell culture | Protein and cellular analyses: |
| - Cell maintenance | - Western blot |
| - Cell counting | - Immunocytochemistry |
| - Cell isolation | - ELISA |

Product information:

	or			
Powder		Liquid	Ambient temperature	48 months

References:

Dulbecco's Phosphate Buffered Saline (DPBS) w/o Calcium w/o Magnesium

 Powder		P0750 - N5L	For 5 liters
		P0750 - N10L	For 10 liters
		P0750 - N50L	For 50 liters
 Liquid	1X	L0615 - 100	100 ml
		L0615 - 500	500 ml
		L0615 - 1000	1000 ml
		L0615 - C10LS	10 liters
 Liquid	10X	X0515 - 100	100 ml
		X0515 - 500	500 ml



For your special formulation, please contact us, we can make your custom-made media on request
Other packagings available on request

Please visit www.biowest.net for more information.

www.biowest.net
Rue de la Caille
49340 Nuaille
Tel.: (+33) 241 464 242
Email: biowest@biowest.net